

Save The Auscultation Sounds Electronic Stethoscope

Model No: HERS® HM-3200



Key Features:





CHARACTERISTIC:

NO MORE TAKING OFF CLOTHES

With a Sound 15 times louder than the traditional Stethoscope, You can hear clear sound even with clothes on, This makes Auscultation for women and the elders easier

Save The Auscultation Sounds

You can save The Auscultation Sounds to your PC and Send information using your E-Mail

Clear The Auscultation Sounds

The Auscultation Sounds can be heard more in detail by minimizing the chafing noise

No Alkaline Battery required for the Stethoscope

We adapted the PC USB Re-Charging System to our Stethoscope

The Power saving mode give the Stethoscope a long lifetime

Enhance your Dignity

The Simple and Up- To- Date design Gives High Satisfaction

SPECIFICATION;

Length: 750mm

Weight; 160g

Battery Lifftime in continue use per charge:

Maximum 38 hours of use.

Amplifies sound up to 15 times.

Two frequency response modes(For Heart sound and Lung sound)

The sounds from the Heart and Lung are maximized while the sounds, from chafing noise is minimized using our circuits

8 level volume control

Save recently used volume level when power off

Power on/off

Micro controller controls automatic power off (Set to two minutes)

Low Battery Indicator

LED Display(Mode, Low Battery. And recharging status)

Speaker provides best sound quality for low frequency sounds

No alkaline battery required

PC USB Interface/Auscultation sounds are automatically stored in DB.

Auscultation sounds can be attached to E-Mail for tele-medicine services

Graphic User Interface based System.

Auscultation sounds are saved and displayed in your PC.

Function Controller:

Power ON/OFF automatic power off after 2 minutes

Volume Control (8 Levels)

Mode Change (Lung mode, Heart Mode)

Display LED

9Mode, Low Battery, Indication, and Recharging status)

Re-Charging using PC USB

Soft Tip(1)/Hard Tip(1)

Human friendly design

Made of silicon which is harmless to the human body

Program CD

Power saving mode embedded

Maximum 38 hours of use

Ear Piece +1

Program CD for program Initialization